

RT-PCR from zebrafish embryos

Elizabeth Wiellette

(also see protocol Grinblat et al, in Methods chapter)

1) Isolate RNA from embryos:

Use about 50 embryos

Use a standard method for isolating RNA, such as TRIZOL (Invitrogen) or RNeasy (Qiagen)

2) DNase treat isolated RNA:

Resuspend the pellet from the isolation step into 100 μ L DNase reaction:

5 μ L (5 U) RNase-free DNase

10 μ L 10X reaction buffer

85 μ L RNase-free H₂O

Incubate at 37°C for 30 minutes

Phenol:CHCl₃ extract and precipitate with NaOAc, glycogen and EtOH

Resuspend pellet in 250 μ L H₂O (not sure about this volume – test possibilities)

3) RT reactions:

Set up two reactions for each RNA sample, for +/- **RT**

1 μ L 25 μ M random hexamer

10 μ L RNA

Incubate at 70°C for 10 min, snap freeze

Add: 4 μ L 5X RT buffer

2 μ L 0.1M DTT

1 μ L 10 mM dNTPs

Incubate at 37°C for 5 minutes

Add: 1 μ L Superscript RT (or do not include this for **-RT** sample)

Incubate at 37°C for 2 hr.

Inactivate at 95°C for 10 minutes

4) PCR amplification of cDNA samples:

Make a common PCR reaction mix with the following for each sample:

37 μL H₂O

5 μL 10X PCR buffer

1 μL 10mM dNTPs

0.5 μL α [³²P]-dCTP (3000 Ci/mmol; 10 mCi/mL; 3.3 μM)

0.5 μL HotStarTaq

Add 44 μL of this mix to:

4 μL template (RT reaction from above)

1 μL each 20 μM primer

Cycles: 15 min. at 95°C

xx cycles of: (cycle number determined empirically for each primer pair)

1 min. at 94°C

1 min. at xx°C (temp determined by primers - probably 55-60°C)

1 min. at 72°C

Note about determining cycle numbers: for RT-PCR to work, the cycle numbers must place the reaction progress within the linear amplification range. This is usually around 20-30 cycles, and to test this, set 4 reactions as above (all with the same primers and template), then remove a tube after 21, 24, 27, 30 cycles.

Primer design: probably want a product that is 300-600 bp

5) Acrylamide gel analysis of PCR products:

Add standard loading buffer to entire 50 μL PCR reaction

Load 17.5 μL of each sample onto a 5% Polyacrylamide gel:

6.25 mL Accugel (40% 19:1 acrylamide:bis)

42.4 mL H₂O

1 mL 50x TAE

0.35 mL 10% APS

17.5 μL TEMED

Run gel at 250V for 1.5-2 hr.

Dry gel at 80°C for 0.5 hr.