

## Single Cell PCR for Zebrafish embryos

1. Prepare PCR tubes. One tube per fish to PCR. Add 5 ul of 10X Qiagen PCR buffer plus 5 ul of ProK/NP40 mix (below) to each tube.
    - a. ProK/NP40 mix
      - i. 1 ul of 20mg/ml ProK
      - ii. 1 ul of 10% NP40
      - iii. 48 ul of H<sub>2</sub>O
  2. Remove small number of cells from blastula stage embryo using micromanipulation/transplantation apparatus in the Bartel Lab.
  3. Put cells into PCR tube with PCR buffer and ProK/NP40.
  4. Incubate tubes in 50 degree water bath for 1 hour.
  5. Make up PCR mix for your desired product.
    - a. PCR mix
      - i. 1ul 10mM dNTP mix
      - ii. Forward Primer of 0.3ul 20uM
      - iii. Reverse Primer of 0.3ul 20uM
      - iv. 2ul 25mM MgCl<sub>2</sub>
      - v. H<sub>2</sub>O up to 50 ul total
      - vi. Taq Polymerase (Qiagen) 0.25ul per reaction  
(don't forget you already have 10 ul in your reaction tube to start with)
- add mix to each PCR reaction tube and start PCR.

Set PCR machine for the desired settings according to your primers and product size.

95 degrees 15 min  
95 degrees 1 min  
60 degrees 1 min (annealing temp should be about 5 degrees lower than primer T<sub>m</sub>)  
72 degrees 1 min per kb product size  
72 degrees 10 min  
4 degrees until you can take it out.

Run on an agarose gel and look for the appropriately sized band.