

Western Blotting

Prepare samples:

For zebrafish embryos:

Option 1:

Take live embryos and put into 1.5 ml tube with E3.

Centrifuge gently for 1-2 minutes

-yolk lipids will rise to the top

-remove supernatant and add sample buffer to embryos and homogenize.

Decide if you need protease/phosphatase inhibitors and add accordingly.

Option 2:

Freeze samples on dry ice/EtOH or in liquid nitrogen.

Thaw and add lysis buffer with protease inhibitors added. (lysis buffer should be cold)

TrisHCl pH 7.4	1ml
EDTA 0.5M	0.4ml
NaCl 5M	3ml
Brij 96 10%	8.75ml
NP40 10%	1.25ml
Bring up to 100 ml with H ₂ O	

Add protease inhibitors and keep on ice after this point

Leupeptin 1:1000

Aprotinin 1:1000

PMSF 1:1000

Add to tissue and homogenize.

After samples are in solution centrifuge in cold for 10 minutes at full speed

Remove the supernatants which contain your sample and keep on ice.

Perform protein assay if desired.

Add 2X SDS sample buffer and proceed.

SDS sample buffer (From Zebrafish book)

0.63 ml 1MTris-HCL, pH6.8

1.0 ml glycerol

0.5 ml Beta-mercaptoethanol

1.75 ml 20% SDS

6.12 ml H₂O

(10 ml total)

add 0..25% bromophenol blue

-store at -20 in aliquots

Pour gel

Set up Biorad (or desired) gel apparatus.

Be sure glass is sealed to prevent leakage.

-pour separating gel (see recipe for desired percentage) using a glass pipet along edge to prevent bubbles from forming in the gel.

-Gently overlay with EtOH and allow to solidify.

-Rinse with H₂O and pour stacking gel.

-Insert desired comb size into stacking gel (be careful of bubbles) and allow to solidify.

-remove combs, set up running apparatus, and fill with running buffer.

Gel Loading/Running

-Take prepared samples including marker proteins

-heat prior to use (heating will prevent/remove precipitates caused by SDS)

-Load samples into well using gel loading tips.

-Set voltage low for the first 20 minutes ~60 Volts, until samples are through the stacking gel layer and lined up evenly at the edge of the separating gel.

-Turn up Voltage to 90-100Volts. Watch sample dye from and stop gel from running once the dye front is at the desired distance through the gel.

Transfer set up

-Before setting up transfer sandwich cut 3MM paper and PVDF to size of gel.

-Wet PVDF in MeOH, rinse with H₂O and soak in transfer buffer.

Soak everything in transfer buffer before putting into sandwich.

Make transfer sandwich as follows

-Sponge

-2 pieces 3MM Paper

-Gel(s)

-PVDF

-2 pieces 3MM Paper

-Sponge

Put sandwich into transfer box, set up voltage box, fill with buffer to completely cover membranes, and run (in the correct direction!!)

Transfer for 1 hour at 100 Volts, or 70-80 for 2 hours, or lower to transfer overnight. (Best to do transfer in the cold room to prevent heating of samples.)

After transfer take apart sandwich and wash PVDF membrane 2X in TBST

-cut corner of PVDF membrane to remember orientation of gel/samples

Blocking and primary antibody incubation

-block in 5% skim milk in TBST for 1 hour at room temperature (or in block desired for best results with your antibody)

-incubate in plastic seal-a-meal bags in primary antibody at desired concentration.

-seal 3 edges,

-add ~3ml of fresh block solution with primary antibody added

-add membrane

-seal without bubbles!

-leave overnight at 4 degrees on shaker

Staining Gel in Coomassie Blue

-carefully put gel into stain for 2 hours on shaker

-destain overnight (see recipes for stain)

Secondary Antibody incubation

After overnight incubation wash 3X in 1XTBST for 30 minutes each

Incubate in the desired concentration of secondary antibody at room temperature in seal-a-meal bag for 2 hours

Wash 3X in 1XTBST for 30 minutes each

Developing Western

-Set up seal-a-meal bag for chemiluminescence reaction with ECL.

-add 1 ml of solution 1.

-add 1 ml of solution 2 and mix well.

-incubate membrane in solution for 1 minute.

-remove and place membrane between two sheet protectors and seal to keep blot moist and flat without bubbles.

-put into film container and add film in dark room.

-expose blot for required amount of time/check initially after 3 minutes and expose accordingly.

-develop in machine on 3rd or 5th floor.

-put blot back in 1XTBST to store in seal-a-meal bag at 4 degrees for stripping and/or reprobing experiments.

Western Recipes

For two THICK 12 % mini-gels

12 ml 30% acrylamide/bis
7.5 ml 4Xtris/SDS pH 8.8
10.5 ml H₂O
100 µl 10% APS
20 µl TEMED

For Stacking gel

1.95 ml 30% acrylamide/bis
3.75 ml 4Xtris/SDS pH 6.8
9.15 ml H₂O
150 µl 10% APS
30 µl TEMED

Block in 0.25% gelatin in TBST

200 ml TBST with 0.5 g gelatin
microwave to dissolve gelatin and cool, should be clear

10X Running buffer

30 g Tris
144 g glycine
10 g SDS
Do not need to pH, should be about 8
Up to 1 liter

Recipe for 10X TBST (for 1 liter)

100 mM Tris Cl pH 7.5	1M Tris	121.1 g
150 mM NaCl	1.5M NaCl	87.66 g
0.1% Tween 20	1% Tween	10 ml

For three liters of 10X TBST

363.3 g Tris Cl

262.98 g NaCl

30 ml Tween 20

up to 3 liters and pH to 7.5 (will take about 100 ml of 10M HCl)

1X BCA lysis buffer with phosphatase inhibitors

25 mM Tris pH 8.0 250 μ l of 1M solution

2 mM EDTA pH 8.0 20 μ l of 0.5 M solution

10% glycerol 1 ml

1% Triton-X 1000 100 μ l

2 mM Na₃VO₄ 100 μ l of 0.2 M solution

20 mM NaF 100 μ l of 2 M solution

Up to 10 ml with H₂O

5X BCA lysis buffer with phosphatase inhibitors

25 mM Tris pH 8.0 1.25 ml of 1M solution

2 mM EDTA pH 8.0 100 μ l of 0.5 M solution

10% glycerol 5 ml

1% Triton-X 1000 500 μ l

2 mM Na₃VO₄ 500 μ l of 0.2 M solution

20 mM NaF 500 μ l of 2 M solution

Up to 10 ml with H₂O

2X Sample Buffer for Protein Gels

25 ml 4X Tris Cl/SDS pH 6.8

20 ml glycerol

2.5 ml β -mercaptoethanol (weigh on good balance)

50 mg bromophenol blue

4 g SDS

Up to 100 ml with H₂O (52.5 ml)

4X TrisCl/SDS

6.05 g Tris base in 40 m (pH 6.8)

Filter 0.45 micron

H₂O to 100 ml

0.4g SDS

4XTrisCl/SDS

91g Tris base 300ml (pH8.8)

H₂O to 500ml total

Filter 0.45 micron

2g SDS

Coomassie Blue Stain

200ml acetic acid

1800 ml H₂O

0.5g Coomassie blue G-250

Mix 1 hr and fliter (No. 1 Whatman paper)

Store at room temp

(final = 0.025% Coomassie blue with 10% acetic acid)

1X Transfer Buffer (4 liters)

12g Tris

57.6g glycine

800ml MeOH

H₂O up to 4 liters